

**CERTIFICATE OF ANALYSIS**  
Swab Products  
**SARS-CoV-2 Positive Swab**

<b>Source Material:</b>	Nasal Swab
<b>Expiry Date:</b>	None assigned. Users should determine suitability in their own system or refer to any known draw date.
<hr/>	
<b>Product Code:</b>	H427-8
<b>Antigen/Marker present:</b>	SARS-CoV-2
<b>Species:</b>	Human
<b>Format:</b>	Dry
<b>Storage Temperature:</b>	≤-15°C. Minimise freeze-thaw cycles
<b>Shipping Temperature:</b>	≤-15°C recommended
<b>Intended Use:</b>	For performance analysis & quality assurance

**Health & Safety:**

---

**Infectious Disease Tests Performed:** The material has not been tested for the presence of infectious agents (other than under “Test Results Data”). See under Precaution for guidance.

**Chemicals:** For further information please request a Safety Data Sheet.

**Precaution:** Handle as hazardous using appropriate Personal Protective Equipment. Refer to Safety Data Sheet (available on request) for safety information and precautions.

**Test Results Data:**

---


**Test performed:** ABI 7500DX CDC Method SARS-CoV-2 PCR (Ct)

**Test results:** See under “Lot/Donor information”.

**Lot/Donor information:**

Lot#	Symptom Date	Collection Date	Donor Age	Donor Gender	PCR Result	N1 (Ct)	N2 (Ct)	RP (Ct)	Swab Type
SW3118-B	28-Jun-2021	01-Jul-2021	32	Male	Positive	31.2215	31.8615	29.1587	Gongdong
SW3074-B	16-Jul-2021	20-Jul-2021	23	Male	Positive	19.3945	21.0624	24.1569	Puritan foam
SW3085-B	02-Jul-2021	05-Jul-2021	53	Male	Positive	26.3652	26.3038	24.2528	Puritan foam
SW3082-B	06-Jul-2021	08-Jul-2021	90	Female	Positive	20.1449	19.4321	28.2528	Puritan foam
SW3084-B	05-Jul-2021	12-Jul-2021	34	Male	Positive	24.5552	25.1585	25.3589	Puritan foam
SW3106-B	03-Jul-2021	08-Jul-2021	22	Female	Positive	34.7407	36.098	24.3554	Puritan foam
SW3113-B	09-Jul-2021	12-Jul-2021	28	Male	Positive	27.7878	27.36	26.713	Puritan foam
SW3111-B	11-Jul-2021	14-Jul-2021	33	Male	Positive	17.4545	16.4514	23.0393	Puritan foam
SW3099-B	13-Jul-2021	14-Jul-2021	25	Female	Positive	18.3992	18.0493	25.4183	Puritan foam
SW3100-B	11-Jul-2021	14-Jul-2021	42	Male	Positive	28.2203	27.6704	26.3609	Puritan foam
SW5131-A	10-Jul-2021	13-Jul-2021	14	Female	Positive	17.3535	17.2477	25.0305	Puritan foam
SW5132-A	17-Jul-2021	20-Jul-2021	35	Male	Positive	17.6671	18.9077	28.0353	Puritan foam
SW2931-B	26-Jul-2021	28-Jul-2021	19	Male	Positive	15.9874	16.112	25.6416	Puritan foam
SW2939-B	16-Jul-2021	22-Jul-2021	14	Male	Positive	21.0676	21.1559	27.0992	Puritan foam
SW2940-B	16-Jul-2021	22-Jul-2021	43	Female	Positive	26.8182	27.2905	24.8832	Puritan foam
SW5133-A	Asymptomatic	27-Jul-2021	27	Male	Positive	37.5112	38.0921	30.1785	Puritan foam
SW5134-A	25-Jul-2021	27-Jul-2021	49	Male	Positive	33.6889	32.5269	32.3963	Puritan foam
SW2957-B	17-Jul-2021	21-Jul-2021	58	Male	Positive	21.1589	22.2602	28.0673	Puritan foam
SW2948-B	11-Jul-2021	13-Jul-2021	11	Female	Positive	25.6192	24.9788	27.2963	Puritan foam
SW5135-A	11-Jul-2021	13-Jul-2021	53	Male	Positive	24.1467	23.7763	32.0131	Puritan foam
SW5136-A	25-Aug-2021	27-Aug-2021	29	Male	Positive	29.3745	31.6026	30.4898	Puritan foam
SW5137-A	30-Aug-2021	31-Aug-2021	45	Male	Positive	35.7064	36.6734	27.2664	Puritan foam
SW3054-B	11-Sep-2021	13-Sep-2021	38	Female	Positive	36.8342	37.4712	22.9407	Puritan foam
SW3056-B	10-Sep-2021	14-Sep-2021	66	Female	Positive	28.3246	30.4086	31.0382	Puritan foam
SW3059-B	12-Sep-2021	14-Sep-2021	38	Female	Positive	21.9476	23.9997	27.8751	Puritan foam
SW3061-B	08-Sep-2021	14-Sep-2021	35	Female	Positive	33.672	36.0364	27.3248	Puritan foam
SW3055-B	10-Sep-2021	13-Sep-2021	19	Female	Positive	33.7483	36.3286	31.2296	Puritan foam
SW3066-B	10-Sep-2021	17-Sep-2021	46	Male	Positive	30.232	31.5866	32.1905	Puritan foam
SW3067-B	15-Sep-2021	20-Sep-2021	21	Male	Positive	36.1805	34.5415	30.2556	Puritan foam
SW3060-B	11-Sep-2021	14-Sep-2021	31	Female	Positive	29.6726	31.5799	32.5708	Puritan foam

**Quality Assurance Approval:**

<b>Signature:</b>	
<b>Date:</b>	17 <sup>th</sup> October 2022

EXAMPLE